

APOBEC3D/APOBEC3F Antibody



PACO00407

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

APOBEC3D/APOBEC3F

Uniprot

Q96AK3/Q8IUX4

Synonyms:

APOBEC3D; Probable DNA dC->dU-editing enzyme APOBEC-3D; APOBEC3F; DNA dC->dU-editing enzyme APOBEC-3F; Apolipoprotein B mRNA-editing enzyme catalytic polypeptide-like 3F

Immunogen:

Synthesized peptide derived from the Internal region of human APOBEC3D/F.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of MCF7 cells using APOBEC3D/F Polyclonal Antibody.

