

AQP4 Antibody



PACO00414

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

AQP4

Uniprot

P55087

Synonyms:

AQP4; Aquaporin-4; AQP-4; Mercurial-insensitive water channel; MIWC; WCH4

Immunogen:

Synthesized peptide derived from the Internal region of human AQP4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of MCF7 cells using AQP4 Polyclonal Antibody.

