

AQP5 Antibody



PACO00415

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

AQP5

Uniprot

P55064

Synonyms:

AQP5; Aquaporin-5; AQP-5

Immunogen:

Synthesized peptide derived from the C-terminal region of human AQP5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HuvEc cells using AQP5 Polyclonal Antibody.

