

ARHGEF10 Antibody



PACO00420

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

ARHGEF10

Uniprot

O15013

Synonyms:

ARHGEF10; KIAA0294; Rho guanine nucleotide exchange factor 10

Immunogen:

Synthesized peptide derived from the Internal region of human ARHGEF10.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of COS7 cells using ARHGEF10 Polyclonal Antibody.

