

ARHGEF19 Antibody



PACO00421

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

ARHGEF19

Uniprot

Q8IW93

Synonyms:

ARHGEF19; Rho guanine nucleotide exchange factor 19; Ephexin-2

Immunogen:

Synthesized peptide derived from the Internal region of human ARHGEF19.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293 cells using ARHGEF19 Polyclonal Antibody.

