

ARRDC2 Antibody



PACO00428

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

ARRDC2

Uniprot

Q8TBH0

Synonyms:

ARRDC2; PP2703; Arrestin domain-containing protein 2

Immunogen:

Synthesized peptide derived from the Internal region of human ARRDC2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293 cells using ARRDC2 Polyclonal Antibody.

