

AVPR1B Antibody



PACO00448

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

AVPR1B

Uniprot

P47901

Synonyms:

AVPR1B; AVPR3; VPR3; Vasopressin V1b receptor; V1bR; AVPR V1b; AVPR V3; Antidiuretic hormone receptor 1b; Vasopressin V3 receptor


Immunogen:

Synthesized peptide derived from the Internal region of human AVP Receptor V3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using AVP Receptor V3 Polyclonal Antibody.