

BARD1 Antibody



PACO00455

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

BARD1

Uniprot

Q99728

Synonyms:

BARD1; BRCA1-associated RING domain protein 1; BARD-1


Immunogen:

Synthesized peptide derived from the N-terminal region of human BARD1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using BARD1 Polyclonal Antibody.