

CNN2 Antibody



PACO00503

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

CNN2

Uniprot

Q99439

Synonyms:

CNN2; Calponin-2; Calponin H2; smooth muscle; Neutral calponin


Immunogen:

Synthesized peptide derived from the Internal region of human Calponin 2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using Calponin 2 Polyclonal Antibody.