

CSE1L Antibody



PACO00514

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF, IP

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, IF:1:200-1:1000,

Protein Background:**Gene ID:**

CSE1L

Uniprot

P55060

Synonyms:

CSE1L; CAS; XPO2; Exportin-2; Exp2; Cellular apoptosis susceptibility protein; Chromosome segregation 1-like protein; Importin-alpha re-exporter


Immunogen:

Synthesized peptide derived from the N-terminal region of human CAS.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using CAS Polyclonal Antibody.