

CEP70 Antibody



PACO00584

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

CEP70

Uniprot

Q8NHQ1

Synonyms:

CEP70; BITE; Centrosomal protein of 70 kDa; Cep70; p10-binding protein


Immunogen:

Synthesized peptide derived from the Internal region of human CEP70.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using CEP70 Polyclonal Antibody.