

CHGA Antibody



PACO00600

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

CHGA

Uniprot

P10645

Synonyms:

CHGA; Chromogranin-A; CgA; Pituitary secretory protein I; SP-I


Immunogen:

Synthesized peptide derived from the Internal region of human Chr-A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Chr-A Polyclonal Antibody.