

CIDEB Antibody



PACO00603

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

CIDEB

Uniprot

Q9UHD4

Synonyms:

CIDEB; Cell death activator CIDE-B; Cell death-inducing DFFA-like effector B


Immunogen:

Synthesized peptide derived from the Internal region of human CIDE-B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAW264.7 cells using CIDE-B Polyclonal Antibody.