

CLCN7 Antibody



PACO00609

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

CLCN7

Uniprot

P51798

Synonyms:

CLCN7; H(+)/Cl(-) exchange transporter 7; Chloride channel 7 alpha subunit; Chloride channel protein 7; CIC-7


Immunogen:

Synthesized peptide derived from the N-terminal region of human CLC-7.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of LOVO cells using CLC-7 Polyclonal Antibody.