

# CRK Antibody



PACO00640

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat, Monkey

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

CRK

**Uniprot**

P46108

**Synonyms:**

CRK; Adapter molecule crk; Proto-oncogene c-Crk; p38

**Immunogen:**


Synthesized peptide derived from human Crk II around the non-phosphorylation site of Y221.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of COS7 cells using Crk II Polyclonal Antibody.