

CCNF Antibody



PACO00654

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

CCNF

Uniprot

P41002

Synonyms:

CCNF; FBX1; FBXO1; Cyclin-F; F-box only protein 1


Immunogen:

Synthesized peptide derived from the C-terminal region of human Cyclin F.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using Cyclin F Polyclonal Antibody.