

CYP4X1 Antibody



PACO00671

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

CYP4X1

Uniprot

Q8N118

Synonyms:

CYP4X1; Cytochrome P450 4X1; CYP4X1


Immunogen:

Synthesized peptide derived from the Internal region of human CYP4X1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using CYP4X1 Polyclonal Antibody.