

TPD52L2 Antibody



PACO00682

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TPD52L2

Uniprot

O43399

Synonyms:

TPD52L2; Tumor protein D54; hD54; Tumor protein D52-like 2


Immunogen:

Synthesized peptide derived from the C-terminal region of human D54.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using D54 Polyclonal Antibody.