

# DMTN Antibody



PACO00698

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:****Gene ID:**

DMTN

**Uniprot**

Q08495

**Synonyms:**

EPB49; DMT; Dematin; Erythrocyte membrane protein band 4.9

**Immunogen:**


Synthesized peptide derived from human Dematin around the non-phosphorylation site of S403.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of 293 cells using Dematin Polyclonal Antibody.