

DLX3 Antibody



PACO00709

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

DLX3

Uniprot

O60479

Synonyms:

DLX3; Homeobox protein DLX-3


Immunogen:

Synthesized peptide derived from the Internal region of human Dlx-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using Dlx-3 Polyclonal Antibody.