

DQX1 Antibody



PACO00724

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

DQX1

Uniprot

Q8TE96

Synonyms:

DQX1; ATP-dependent RNA helicase DQX1; DEAQ box polypeptide 1


Immunogen:

Synthesized peptide derived from the C-terminal region of human DQX1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using DQX1 Polyclonal Antibody.