

DTYMK Antibody



PACO00726

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

DTYMK

Uniprot

P23919

Synonyms:

DTYMK; CDC8; TMPK; TYMK; Thymidylate kinase; dTMP kinase


Immunogen:

Synthesized peptide derived from the N-terminal region of human dTMP Kinase.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using dTMP Kinase Polyclonal Antibody.