

DNAL1 Antibody



PACO00729

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

DNAL1

Uniprot

Q4LDG9

Synonyms:

DNAL1; C14orf168; Dynein light chain 1; axonemal


Immunogen:

Synthesized peptide derived from the Internal region of human Dynein LC 1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MCF7 cells using Dynein LC 1 Polyclonal Antibody.