

E2F4/E2F5 Antibody



PACO00732

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

E2F4/E2F5

Uniprot

Q16254/Q15329

Synonyms:

E2F4; Transcription factor E2F4; E2F-4; E2F5; Transcription factor E2F5; E2F-5


Immunogen:

Synthesized peptide derived from the Internal region of human E2F-4/5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using E2F-4/5 Polyclonal Antibody.