

ELOVL4 Antibody



PACO00756

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ELOVL4

Uniprot

Q9GZR5

Synonyms:

ELOVL4; Elongation of very long chain fatty acids protein 4; 3-keto acyl-CoA synthase
ELOVL4; ELOVL fatty acid elongase 4; ELOVL FA elongase 4


Immunogen:

Synthesized peptide derived from the Internal region of human ELOVL4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using ELOVL4 Polyclonal Antibody.