

RASA3 Antibody



PACO00836

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

RASA3

Uniprot

Q14644

Synonyms:

RASA3; Ras GTPase-activating protein 3; GAP1(IP4BP); Ins P4-binding protein


Immunogen:

Synthesized peptide derived from the Internal region of human GAP1-InsP4 BP.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using GAP1-InsP4 BP Polyclonal Antibody.