

GAP43 Antibody



PACO00837

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

GAP43

Uniprot

P17677

Synonyms:

GAP43; Neuromodulin; Axonal membrane protein GAP-43; Growth-associated protein 43; Neural phosphoprotein B-50; pp46


Immunogen:

Synthesized peptide derived from human GAP-43 around the non-phosphorylation site of S41.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HELA KB cells using GAP-43 Polyclonal Antibody.