

GATA4 Antibody



PACO00841

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

GATA4

Uniprot

P43694

Synonyms:

GATA4; Transcription factor GATA-4; GATA-binding factor 4


Immunogen:

Synthesized peptide derived from human GATA-4 around the non-phosphorylation site of S105.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using GATA-4 Polyclonal Antibody.