
Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

GK

Uniprot

P32189

Synonyms:

GK; Glycerol kinase; GK; Glycerokinase; ATP: glycerol 3-phosphotransferase


Immunogen:

Synthesized peptide derived from the C-terminal region of human GK1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using GK1 Polyclonal Antibody.