

GPR176 Antibody



PACO00872

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

GPR176

Uniprot

Q14439

Synonyms:

GPR176; Probable G-protein coupled receptor 176; HB-954


Immunogen:

Synthesized peptide derived from the C-terminal region of human GPR176.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using GPR176 Polyclonal Antibody.