

GSTT1/GSTT4 Antibody



PACO00892

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

GSTT1

Uniprot

P30711/A8MPT4

Synonyms:

GSTT1; Glutathione S-transferase theta-1; GST class-theta-1; Glutathione transferase T1-1; GSTT4; Glutathione S-transferase theta-4; GST class-theta-4


Immunogen:

Synthesized peptide derived from the N-terminal region of human GSTT1/4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using GSTT1/4 Polyclonal Antibody.