

HBP1 Antibody



PACO00896

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

HBP1

Uniprot

O60381

Synonyms:

HBP1; HMG box-containing protein 1; HMG box transcription factor 1; High mobility group box transcription factor 1


Immunogen:

Synthesized peptide derived from the C-terminal region of human HBP1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using HBP1 Polyclonal Antibody.