

HIPK4 Antibody



PACO00906

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

HIPK4

Uniprot

Q8NE63

Synonyms:

HIPK4; Homeodomain-interacting protein kinase 4


Immunogen:

Synthesized peptide derived from the C-terminal region of human HIPK4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using HIPK4 Polyclonal Antibody.