

HLA-DOB Antibody



PACO00916

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

HLA-DOB

Uniprot

P13765

Synonyms:

HLA-DOB; HLA class II histocompatibility antigen; DO beta chain; MHC class II antigen
DOB


Immunogen:

Synthesized peptide derived from the N-terminal region of human HLA-DOBeta.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using HLA-DOBeta Polyclonal Antibody.