

# HOXB5 Antibody



PACO00932

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

HOXB5

**Uniprot**

P09067

**Synonyms:**

HOXB5; HOX2A; Homeobox protein Hox-B5; Homeobox protein HHO. C10; Homeobox protein Hox-2A; Homeobox protein Hu-1

**Immunogen:**


Synthesized peptide derived from the Internal region of human HoxB5.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of HuvEc cells using HoxB5 Polyclonal Antibody.