

IFNG Antibody



PACO00948

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

IFNG

Uniprot

P01579

Synonyms:

IFNG; Interferon gamma; IFN-gamma; Immune interferon


Immunogen:

Synthesized peptide derived from the Internal region of human IFN-gamma.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using IFN-gamma Polyclonal Antibody.