

IKZF3 Antibody



PACO00951

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

IKZF3

Uniprot

Q9UKT9

Synonyms:

IKZF3; ZNFN1A3; Zinc finger protein Aiolos; Ikaros family zinc finger protein 3


Immunogen:

Synthesized peptide derived from the Internal region of human Ikaros 3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using Ikaros 3 Polyclonal Antibody.