

INSL4 Antibody



PACO00968

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

INSL4

Uniprot

Q14641

Synonyms:

INSL4; Early placenta insulin-like peptide; EPIL; Insulin-like peptide 4; Placentin


Immunogen:

Synthesized peptide derived from the Internal region of human INSL4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HT29 cells using INSL4 Polyclonal Antibody.