

JAK2 Antibody



PACO00988

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

JAK2

Uniprot

O60674

Synonyms:

JAK2; Tyrosine-protein kinase JAK2; Janus kinase 2; JAK-2


Immunogen:

Synthesized peptide derived from human JAK2 around the non-phosphorylation site of Y570.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using JAK2 Polyclonal Antibody.