

# JMY Antibody



PACO00989

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

JMY

**Uniprot**

Q8N9B5

**Synonyms:**

JMY; Junction-mediating and -regulatory protein

**Immunogen:**


Synthesized peptide derived from the C-terminal region of human JMY.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of COLO205 cells using JMY Polyclonal Antibody.