

TOR1AIP1 Antibody



PACO01013

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TOR1AIP1

Uniprot

Q5JTV8

Synonyms:

TOR1AIP1; Torsin-1A-interacting protein 1; Lamin-associated protein 1B; LAP1B


Immunogen:

Synthesized peptide derived from the Internal region of human LAP1B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using LAP1B Polyclonal Antibody.