

PACO01014

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

CERS4

**Uniprot**

Q9HA82

**Synonyms:**

CERS4; LASS4; Ceramide synthase 4; CerS4; LAG1 longevity assurance homolog 4

**Immunogen:**


Synthesized peptide derived from the Internal region of human LASS4.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of Jurkat cells using LASS4 Polyclonal Antibody.