

PACO01017

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

LATS1/LATS2

Uniprot

O95835/Q9NRM7

Synonyms:

LATS1; WARTS; Serine/threonine-protein kinase LATS1; Large tumor suppressor homolog 1; WARTS protein kinase; h-warts; LATS2; KPM; Serine/threonine-protein kinase LATS2; Kinase phosphorylated during mitosis protein; Large tumor suppressor ho

Immunogen:

Synthesized peptide derived from human LATS1/2 around the non-phosphorylation site of T1079/1041.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of RAT-MUSCLE cells using LATS1/2 Polyclonal Antibody.

