

MED18 Antibody



PACO01054

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

MED18

Uniprot

Q9BUE0

Synonyms:

MED18; Mediator of RNA polymerase II transcription subunit 18; Mediator complex subunit 18; p28b

Immunogen:

Synthesized peptide derived from the Internal region of human Med18.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of Jurkat cells using Med18 Polyclonal Antibody.

