

MN1 Antibody



PACO01075

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

MN1

Uniprot

Q10571

Synonyms:

MN1; Probable tumor suppressor protein MN1


Immunogen:

Synthesized peptide derived from the Internal region of human MN1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using MN1 Polyclonal Antibody.