

MRPL15 Antibody



PACO01092

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

MRPL15

Uniprot

Q9P015

Synonyms:

MRPL15; HSPC145; 39S ribosomal protein L15; mitochondrial; L15mt; MRP-L15


Immunogen:

Synthesized peptide derived from the Internal region of human MRP-L15.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MCF7 cells using MRP-L15 Polyclonal Antibody.