

MRPL20 Antibody



PACO01095

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

MRPL20

Uniprot

Q9BYC9

Synonyms:

MRPL20; 39S ribosomal protein L20; mitochondrial; L20mt; MRP-L20


Immunogen:

Synthesized peptide derived from the C-terminal region of human MRP-L20.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using MRP-L20 Polyclonal Antibody.