

MRPL21 Antibody



PACO01096

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

MRPL21

Uniprot

Q7Z2W9

Synonyms:

MRPL21; 39S ribosomal protein L21; mitochondrial; L21mt; MRP-L21


Immunogen:

Synthesized peptide derived from the Internal region of human MRP-L21.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using MRP-L21 Polyclonal Antibody.