

MTHFR Antibody



PACO01117

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

MTHFR

Uniprot

P42898

Synonyms:

MTHFR; Methylenetetrahydrofolate reductase


Immunogen:

Synthesized peptide derived from the Internal region of human MTHFR.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using MTHFR Polyclonal Antibody.