

MVK Antibody



PACO01121

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

MVK

Uniprot

Q03426

Synonyms:

MVK; Mevalonate kinase; MK


Immunogen:

Synthesized peptide derived from the Internal region of human MVK.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using MVK Polyclonal Antibody.